

Daily patient setup made easy

Position of the cushions is easily reproduced during each fraction. Comprehensive color coding facilitates the use of the cushion combinations.

Make patient comfort a priority

Minimal pressure points and soft, ergonomically designed cushions.

Made to last

A busy radiation therapy department requires high durability devices to withstand rough handling conditions.



Maximize clearance, minimize collision

The full AIO setup fits in small CT-bores and reduces the risk of collision with the gantry.

Immobilize without a base plate


Cushions can be positioned directly on the treatment table without the use of a base plate.



The AIO Solution®

With one base plate, four cushion sets and thermoplastic masks, it is possible to position and immobilize patients with lesions in brain, head and neck, thorax, lung and pelvis, in a precise, reproducible and comfortable way.

Prone breastboard and belly- and pelvicboard cushion sets are under development.



About ORFIT

Orfit brings high precision and comfort to the positioning and immobilization of cancer patients in Radiation Oncology. The combination of Orfit thermoplastic masks and positioning devices provides optimal stability to a point where daily reproducibility and patient movement is limited to less than 1.5 millimeter.

Discover the complete family of Orfit systems for brain, head and neck, supine and prone breast treatment, pelvis and abdomen treatment, SRS, SBRT, extremities, pediatric, MR and proton therapy on www.orfit.com



The AIO Solution®

Visit
WWW.ORFIT.COM/AIO
for more information.

www.orfit.com



ORFIT INDUSTRIES
Vosveld 9A | B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

03-04-2019
Art. N° 51014/1

orfit The precision expert